

TEM FIXATION PROTOCOLS:

Chemicals needed

Post-Fixative:

1% Osmium Tetroxide:

- 0.1g OsO₄ Crystals
- 10 mL DH₂O
- or
- 2 mL 4% OsO₄
- 4 mL Sorensen's Buffer
- 2 mL DH₂O (can substitute buffer if desired)

Dehydration Series:

Ethanol: 10%, 30%, 50%, 70%, 90%, 100%, 100% in Zeolite beads (molecular sieve)

Acetone: 10%, 30%, 50%, 70%, 90%, 100%, 100% in Zeolite beads (molecular sieve)

Protocol (derivations exist depending on your sample):

1. Under a hood and while wearing gloves, place fixative of choice

9. Prepare the embedding resin according to manufacturer instructions. Make small aliquots of the following ratios:

Spurr's:

1. Spurr's/acetone: 25/75

2. Spurr's/acetone: 50/50

3. Spurr's/acetone: 75/25

4. Spurr's – 100%

London Resin:

LR/ethanol: 25/75

LR/ethanol: 50/50

LR/ethanol: 75/25

LR-